cerillo
A HAMILTON Company



Canopy Technical Specifications

The Canopy is RoHS, REACH and FCC compliant, CE certified, made in the USA, and has a two-year warranty.

General

Compatible Products	Cerillo Microplate Readers
Dimensions	11.32 X 6.61 X 4.78 cm (LWH)
Weight	170 g
Power Input	5.1 V, 1 A via USB-C receptacle

Communication

Wireless Protocol	IEEE 802.15.4
Wireless Range	15 m

Cerillo Labrador Premium Software

Operating System	macOS: 12 (Monterey); Windows: 10, 11

All Cerillo products are for research use only and are not intended for use in diagnostic procedures.
Contact us: www.cerillo.bio | info@cerillo.bio